




Alexis[®]

Orthopaedic Protector



Alexis


Orthopaedic Protector

Shields soft tissue – including skin, fat, muscles and nerves – from contact with abrasive instrumentation and debris at the wound margins.



Distributes the force of retraction evenly for optimal exposure of the joint capsule.

Provides 360 degrees of atraumatic protection and retraction.

Total Hip Arthroplasty



Total Shoulder Arthroplasty





Alexis Orthopaedic Protector

Featuring a rigid retraction ring for maximum retraction and a flexible retraction ring for maximum versatility

Reorder No.	Size	Sheath Length	Incision Range	Qty/Box
-------------	------	---------------	----------------	---------

Rigid Retraction Ring

HR000	Extra Small/Medium	14cm	2.5–7cm	5
HR001	Small/Small	14cm	2.5–8cm	5
HR004	Small/Medium	14cm	2.5–8cm	5
HR005	Medium/Large	17cm	5–13cm	5

Flexible Retraction Ring

HR100	Extra Small/Medium	14cm	2.5–7cm	5
HR101	Small/Small	14cm	2.5–8cm	5
HR104	Small/Medium	14cm	2.5–8cm	5
HR105	Medium/Large	17cm	5–13cm	5

Visit [appliedmedical.com/alexisortho](https://www.appliedmedical.com/alexisortho) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2023 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

11/2023 320451-EN-C

