



Epix™

Laparoscopic Clip Appliers

Epix

Universal Clip Applier

A Universal Solution for Ligation

The Epix universal clip applier delivers a true M/L clip through a 5mm shaft when the device is used along with a Kii™ 5mm trocar.



True M/L clip


TiGold Clip Technology

Contributes to better clip retention by subjecting the clips to a proprietary heat treatment process.



Low-Profile Jaw

Facilitates visualization and provides clearance around the tissue to be ligated.



Note: Numbers are based on an internal test report from May 2023.



Epix Laparoscopic Clip Appliers

Reorder No.	Description	Shaft Diameter	Qty/Box
CA500 ^a	Epix universal clip applier	5mm	3
CA700 ^b	Epix 10mm clip applier	10mm	3

Note: All Epix clip appliers deliver 20 M/L TiGold clips.

^a Model CA500 is compatible with Applied Medical Kii 5mm trocars and larger trocars.

^b Model CA700 is compatible with all 10mm and larger trocars.

Devices listed must be administered by a healthcare professional. Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2025 Applied Medical Resources Corporation, Rancho Santa Margarita, California. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ™ are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

All other product names, logos, brands, and registered or unregistered trademarks referenced herein are the property of their respective owners, and no trademark rights to the same are claimed. Use of these names, trademarks and brands does not imply endorsement.

07/2025 666393-EN-A

