



Inzii[®]
Retrieval Systems

SPECIMEN RETRIEVAL BAG OFFERINGS

Inzii retrieval systems, with bags available in either polyurethane material or ripstop nylon material, are designed to facilitate specimen extraction through small linear defects. The innovative guide bead makes it easy to reopen and close the bag, while the extra-large opening allows for smooth delivery of the specimen into the bag. Offering clinical versatility, durability and ease of use, Inzii retrieval systems maximize specimen security during tissue removal.



**Inzii Retrieval Systems
In-Service**



**Inzii
Ripstop Redeployable
Retrieval System
In-Service**


Inzii

Retrieval Systems (Polyurethane Bags)



Inzii 12/15mm

15mm bag in a 12mm device



Inzii 10mm

Standard specimen retrieval bag



Inzii universal

10mm bag in a 5mm* device

* Applied Medical Kii trocar (if using a non-Kii trocar, a larger port size will be needed)

Inzii

Retrieval Systems (Polyurethane Bags)



Guide Bead and Cord Loop

Allow for reopening and closing the bag and facilitate specimen extraction through small linear defects.



Deployment System


Maximizes clinical versatility and minimizes the need for large ports.



Reinforced Bag

Ensures specimen security and minimizes the risk of contamination.

Burst Strength Comparison



Note: Data is on file.

Inzii

Ripstop Redeployable Retrieval System




Guide Bead and Cord Loop

Allow for reopening and closing the bag after final deployment and facilitate specimen extraction through small linear defects.



Redeployable System

Enables multiple deployments and multiple specimen extractions within the same procedure.



Ripstop Nylon Bag Material

Provides strength and aids in secure specimen retrieval.



Inzii Universal Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD003	5mm ^a	180mL	10

^a Applied Medical Kii trocar (if using a non-Kii trocar, a larger port size will be needed)



Inzii 10mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD001	10mm	225mL	10



Inzii Ripstop Redeployable Retrieval System

Reorder No.	Port Size	Volume	Qty/box
CD005	11mm	265mL	5



Inzii 12/15mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD004	12mm	1,600mL	5

Visit [appliedmedical.com/inzii](https://www.appliedmedical.com/inzii) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2024 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union. All other product names, logos, brands, registered or unregistered trademarks referenced herein are the property of their respective owners and no trademark rights to the same are claimed. Use of these names, trademarks and brands does not imply endorsement.

01/2024 411131-EN-C

