



Applied Scope Warmer

The background of the slide is white with a large teal geometric shape on the left side. This shape is a large triangle that points towards the bottom right corner, with a thin, slightly darker teal line running parallel to its hypotenuse.

A SIMPLE SOLUTION FOR BETTER VISIBILITY

The Answer to Laparoscopic Fogging

The Applied scope warmer retains heat throughout lengthy procedures, eliminates the need for chemical gels or wipes, and minimizes the need for repeated scope cleansing and defogging.

Scope-Warmer Seal

Retains heat throughout lengthy procedures.

Scope-Warmer Bottle

Minimizes fogging by prewarming the laparoscope to body temperature.

Scope-Warmer Base

Stabilizes the system when a scope is inserted.





Applied Scope Warmer

Reorder No.	Description	Qty/Box
C3001	Scope-warmer bottle ^a Reusable Volume capacity = 500mL	1
C3002	Scope-warmer base ^a Reusable	1
C3101	Scope-warmer seals Disposable, single use Individually packaged in sterile pouches	20

Note: The scope is not included.

^a Recommended sterilization method: Autoclave or STERRAD 100S system.

The Applied scope warmer prevents fogging by pre-warming the laparoscope to body temperature.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2025 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ™ are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

All other product names, logos, brands, registered or unregistered trademarks referenced herein are the property of their respective owners and no trademark rights to the same are claimed. Use of these names, trademarks and brands does not imply endorsement.

01/2025 411033-EN-B

