



Applied Insufflation Needle

A large teal graphic element consisting of two overlapping triangles. One triangle is a darker shade of teal and points towards the bottom right corner. The other is a lighter shade of teal and points towards the top left corner, creating a diagonal line across the page.

ENABLING HIGH-FLOW INSUFFLATION

13-gauge percutaneous needle

Tactile and visual indicators

Ergonomic housing design

Applied Insufflation Needle



Spring-Loaded, Blunt Stylet with Red Indicator Button

Provides tactile and visual feedback when the stylet is retracted.

13-Gauge High-Flow Needle

Offers faster insufflation than traditional needles, reducing operating room time.



Ergonomic Housing

Ensures control during needle insertion.





Applied Insufflation Needle

Reorder No.	Size	Qty/Box
C2201	13 gauge, 120mm	20
C2202	13 gauge, 150mm	20

Visit [appliedmedical.com/laparoscopic](https://www.appliedmedical.com/laparoscopic) for more information.

The Applied insufflation needle is used for establishing pneumoperitoneum prior to laparoscopic surgery.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2025 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ™ are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

03/2025 411136-EN-C

