




# Stealth<sup>®</sup>

Atraumatic Occlusion System

## Advantages of Stealth Clips

Stealth clips are designed to provide occlusion where space is at a premium.



### Reduced Spring Force

Three options allow for 1/4-force, 1/2-force or 3/4-force occlusion, depending on the surgeon's preference and the needs of each procedure.



### *latis*<sup>®</sup>/*latis* Clips

Reorder No.	Description	Size
A1601	<i>latis/latis</i> clip	6mm (1/2-force)
A1602	<i>latis/latis</i> clip	6mm (1/4-force)
A1603	<i>latis/latis</i> clip	6mm (3/4-force)



### Soft/Soft Clips

Reorder No.	Description	Size
A1604	Soft/Soft clip	6mm (3/4-force)



## Soft/Fibra® Clips

Reorder No.	Description	Size
A1652	Soft/Fibra clip	6mm (1/2-force)
A1654	Soft/Fibra clip	6mm (1/4-force)
A1650	Soft/Fibra clip	6mm (3/4-force)




## Stealth Surgical Clip Applicator

Reorder No.	Description	Handle Angle	Overall Length
A3200	Stealth surgical clip applicator	13°	7.25"

## Advantages of the Stealth Occlusion System

The low-profile Stealth system has become a leader in atraumatic occlusion because of its superior traction, reduced occlusion force and minimized clip migration.

## Stealth Clamps Versus Standard Clamps



## Ease of Installation

The inserts attach securely and easily.



## Variety of Sizes

Stealth systems offer a wide variety of clamp insert sizes and configurations.

Size 1, 30mm



Size 2, 45mm




Size 3, 60mm



Size 5, 90mm



Sidebite, small  
50mm



Sidebite, medium  
75mm

*Actual sizes*

## Stealth Inserts

Stealth clamp inserts (all of them non-latex) are available in several configurations.

**Soft** and **Traction** inserts are paired with **latis** and **Fibra** inserts to meet specific occlusion requirements.

Clinically proven **latis** pads gently engage the tissue, reducing the pressure necessary to maintain traction. The **Fibra** insert is covered with a woven layer of flexible, soft, finger-like fibers. These fibers mesh with the adventitia to provide traction without added pressure.



### Soft/Fibra Inserts

Reorder No.

Size

*Straight*

A0G01 30mm

A0G02 45mm

A0G03 60mm

A0G04 90mm



## Soft/Traction Inserts

Reorder No. Size

*Straight*

A0A01	30mm
A0A02	45mm
A0A03	60mm
A0A04	90mm



## Soft/Soft Inserts

Reorder No. Size

*Straight*

A0A11	30mm
A0A12	45mm
A0A13	60mm
A0A14	90mm



## Fibra/Fibra Inserts

Reorder No. Size

*Straight*

A0G11	30mm
A0G12	45mm
A0G13	60mm
A0G14	90mm

*Sidebites*

A0G05	Small
A0G06	Medium



## latis/latis Inserts

Reorder No. Size

*Straight*

A0C01	30mm
A0C02	45mm
A0C03	60mm
A0C04	90mm

*Sidebites*

A0C05	Small
A0C06	Medium

## Stealth Clamps

Atraumatic Stealth clamps offer an unmatched level of versatility.

### Consistent, Gentle Occlusion


The jaws provide consistent, gentle occlusion as a result of their parallel alignment.

### Low Profile

The low-profile design enables optimal occlusion in hard-to-access areas.

### Firmly Secured Inserts

Innovative attachment technology fastens the insert securely to the clamp.





### Ensured Durability

Stainless steel construction ensures durability.

### Enhanced Stability

The exceptionally crafted, fine-tuned box lock enhances stability.

## Preferred Configurations



### Bahnsen

The Stealth system's best-selling model



### Wylie

Deep reach with a unique bend



### Sidebite

The market's only padded sidebite

## Stealth Clamp Configurations



Reorder No.	Size	Handle Angle	Overall Length	Jaw Angle	Jaw Length	Applications
A3201	1	15°	6.5"	Straight	30mm	Femoral and iliac
A3202	1	25°	6"	60°	30mm	Femoral and iliac

## Stealth Clamp Configurations



Reorder No.	Size	Handle Angle	Overall Length	Jaw Angle	Jaw Length	Applications
A3203	2	15°	7"	Straight	45mm	Femoral and iliac
A3204	2	15°	7"	34°	45mm	Femoral and iliac
A3205	2	Straight	6.5"	65°	45mm	Femoral and iliac
A3206	2	15°	6.5"	65°	45mm	Femoral and iliac
A3207	3	10°	11.5"	Straight	60mm	Iliac and abdominal aortic
A3208	3	15°	10"	60°	60mm	Iliac and abdominal aortic
A3209	3	35°	10"	60°	60mm	Iliac and abdominal aortic
A3210	3	35°	10"	Bahnson	60mm	Iliac and abdominal aortic
A3211	3	60°	10.5"	Wylie	60mm	Iliac, abdominal aortic, femoral and renal
A3212	5	15°	12"	Straight	90mm	Aortic cross-clamp
A3213	5	35°	11"	60°	90mm	Aortic cross-clamp
A3214	5	35°	10.5"	Bahnson	90mm	Aortic cross-clamp
A3215	Small	25°	7"	Sidebite	50mm	Aortic and vena cava
A3216	Medium	45°	8"	Sidebite	75mm	Aortic and vena cava
A3219	3	35°	8.5"	Bahnson	60mm	Aortic cross-clamp
A3220	5	35°	9"	Bahnson	90mm	Aortic cross-clamp
A3221	2	8.5°	7.5"	Debakey	45mm	Femoral and iliac

Visit [appliedmedical.com/stealth](https://www.appliedmedical.com/stealth) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

08/2022 213951-EN-D

